



Sterilucent™ Sterilization Container System



The Sterilucent Sterilization Container System is built to provide long-lasting protection and sterility maintenance for your valuable instrument investment. Sterilucent containers are made of anodized aluminum and are ergonomically designed with large handles to encourage proper lifting.

Sterilucent sterilization containers are available in multiple sizes and can be configured with a large variety of accessories for use with the Sterilucent PSD-85 Hydrogen Peroxide Sterilizer.

Sterilucent, Inc.
1400 Marshall Street NE
Minneapolis, MN 55413
+1 612 767 3260
+1 877 721-8405
www.sterilucent.com